

Act Fast™

Anti Choking Trainer



Your Anti Choking Trainer should include:


1. The patented Act Fast Vest.
2. An internal Air Bladder
3. Clear Plastic Airway (this should fit snugly into the coupler on the air bladder).
4. Three packages of reusable Foam Plugs, 1 green and 2 yellow.
5. If you have the Red Back Slap Pad model, this should include both a plastic plate and foam pad inside.
6. User Manual, please refer to for detailed instructions.
7. Nylon Drawstring Carry Bag.

IMPORTANT:

First try the Anti Choking Trainer on yourself to ensure it is working properly before using on another person.

It should not require force to use. The air bladder should be soft enough that proper technique will easily expel the foam plug.

For any questions please contact Act Fast.



FOAM PLUGS:

To customize the level of ease for your students' needs, the foam plugs can be rolled in your fingertips to soften, or trimmed to shorten. Either method will allow the foam plug to be expelled more easily if needed.

Limited Warranty

Act Fast warrants to the Customer that the Equipment, (constituting the Anti Choking Trainer) will be free from defects in material and workmanship under normal use and service for a period of one (1) year from date of shipment.

During such period Act Fast will at no charge to the Customer, either repair or replace (at Act Fast's sole option) any part of the Equipment found to be defective in material or workmanship. Act Fast will pay all shipping costs.

Act Fast shall not be responsible for any Equipment defect, the failure of the Equipment to perform any specified function,

or any other nonconformance of the Equipment caused by or attributable to: (i) any modification of the Equipment by the Customer, (ii) the use of the Equipment with any associated or complementary Equipment, or accessory not supplied by Act Fast; (iii) any misuse or abuse of the Equipment; or (iv) installation or use of the Equipment other than in accordance with Act Fast's instructions.

info@actfastmed.com
warranty@actfastmed.com
or
Call toll free 855-934-9340